



# Recombinant Human Flt3-interacting zinc finger protein 1 (FIZ1)

<b>Product Code</b>	CSB-EP839415HU-B
<b>Abbreviation</b>	FIZ1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96SL8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDVPAPTPA PAPPAAAAPR VPFHCSECGK SFRYRSDLRR HFARHTALKP HACPRCGKGF KHSFNLANHL RSHTGERPYR CSACPKGFRD STGLLHHQVV HTGEKPYCCL VCELRFSRS SLGRHLKRQH RGVLPSP LQP GPGLPALSAP CSVCCNVGPC SVCGGSGAGG GEGPEGAGAG LGSWGLAEAA AAAAASLPPF ACGACARRFD HGRELAHWA AHTDVKPFKC PRCERDFNAP ALLERHKLTH DLQGP GAPP QAWAAGPGAG PETAGEGTAA EAGDAPLASD RRLLLGPAGG GVPKLG LLLP EGGGEAPAPA AAAEPS EDTL YQCDCG TFFA SAAALASHLE AHSGPATYGC GHCGALYAAL AALEEHRVS HEGGGGEEAA TAAREREPAS GEPPSGSGRG KKIFGCSECE KLFRSPRDLE RHVLVHTGEK PFPCLECGKF FRHECYLKRH RLLHGTERPF PCHICGKGFI TLSNLSRHLK LHRGMD
<b>Source</b>	E.coli
<b>Target Names</b>	FIZ1
<b>Protein Names</b>	Recommended name: Flt3-interacting zinc finger protein 1 Alternative name(s): Zinc finger protein 798
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes zinc finger protein, which interacts with a receptor tyrosine kinase involved in the regulation of hematopoietic and lymphoid cells. This gene product also interacts with a transcription factor that regulates the expression of rod-specific genes in retina.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.