



Recombinant Schizosaccharomyces pombe Meiotically up-regulated gene 93 protein (mug93)

Product Code	CSB-MP839432SXV
Abbreviation	mug93
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96WW2
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MEEKENVNLV EKSNYVALEN TREIDVFDEF LNAIGNENTI TPVYADSSLT HLRKKSYTKV VHDCTYALVI NPYDKKVIWR RGLAYLRLGH PHLANRDWEH SLELDPNNTY IQKSLHRLKE VYIYRECAE TWQLRHLRVA SSQQLPVGLR KQYPNIIRGK IWWKKVHDNC QLCGQQCELE KENLSAMRSM LYMANTYAKD DTENHSPSAQ IGIESSEDEL ENKITKGEHS LLVPEELYRS NYPCPQNIDQ FLYMIKVLSA PCLYIETFSF PISTINQLFK AHGMSVEQLN LFLKSIHYIG LCSRFCQWS DKARSLMQAL SGLPWFSFVV QHCLHITAAQ ILLHIPDIQE EEFRNWHVSK KPINNTDLSS EFEIAEIPIN CYT
Source	Mammalian cell
Target Names	mug93
Protein Names	Recommended name: Meiotically up-regulated gene 93 protein
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.