



# Recombinant *Thermoplasma volcanium* Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-EP839488TDL
<b>Abbreviation</b>	srp54
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q979Y8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Thermoplasma volcanium</i> (strain ATCC 51530 / DSM 4299 / JCM 9571 / NBRC 15438 / GSS1)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVLESLSSSL RETIRKITGS TYIDKNTVKE ISKDLQRILL KADVNVQTVL NITKEMERRA LEEKPPVGMA HQDYMVKIY EELLKILGEP SNVELKPQTI MLVGLYGHGK TTSAGKLSRF FMKKGLSAGL IAADVHRYAA YEQLQQIAES VNAKFYGDPAK EKDPVKIKN GLEVLKDAQV KIIDTSGRDS LDSDLIEEIR RIKEAVKPDQ VLMVIDATMG QQAGPESKAI NDAVGVTGII ITKMDGTAKG GGALSAVAKM RVPIYFIGTG EHMDDMEIFD PPKFLSRLLG LGDLETLFET VKEANITEEE AQKSFELMT GKFNLKDMYD VWEKFSRPLG MKKLVDSLPL AKLPGSQKIN DDQILTAEQK LQRYRIIMDS MTFQELENPE IINAKRITRI ARGAGVREED VRMLLKEFKA MKNNMKLMKG NRGFKKILQS NFRSGNFGLD DLGLK
<b>Source</b>	E.coli
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-455
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.