



Recombinant Arabidopsis thaliana E2F transcription factor-like E2FE (E2FE)

Product Code	CSB-EP839710DOA-B
Abbreviation	E2FE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LSZ4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSDLSPERFK LAVTSPSSIP ESSLALQLHH SYSRKQKSLG LLCTNFLALY NREGIEMVGL DDAASKLGVE RRRIYDIVNV LESVGLVTRR AKNQYTWKGF SAIPGALKEL QEEGVKDTFH RFYVNEVKG SDEDDDEES SQPHSSSQTD SSKPGSLPQS SDPSKIDNRR EKSLGLLTQN FIKLFICSEA IRIISLDDAA KLLLGDAHNT SIMRTKVRRL YDIANVLSSM NLIKHTHTLD SRKPAFKWLG YNGEPTFTLS SDLLQLESRK RAFGTDITNV NVKRSKSSSS SQENATERRL KMKKHSTPES SYNKSFDVHE SRHGSRGGYH FGPFPAGTGT YPTAGLEDNS RRAFDVENLD SDYRPSYQNQ VLKDLFSHYM DAWKTWFSEV TQENPLPNTS QHR
Source	E.coli
Target Names	E2FE
Protein Names	Recommended name: E2F transcription factor-like E2FE Alternative name(s): DP-E2F-like protein 1 E2F-like repressor E2L3
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.