



Recombinant *Drosophila simulans* Chitinase-like protein Idgf1 (Idgf1)

Product Code	CSB-BP839766DMJ
Abbreviation	Idgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8MX41
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila simulans</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	ANNLVCYYDS TSYLRQGLAK MHTNELDLAL QFCTHLVYGY AGLKSGTLEL FSLNVDLDMF YYKDITALRQ KFPQLKILLS VGGDRDVDEA HPNKYVELLE ANRTAQQNF I DSSMILLKRN GFDGLDLAFQ LPRNKPRKVH GSLGTYWKSF KKLFTGDFVV DPLAEQHKSQ FTDLVGNIKN AFRSANLMLS LTVLPNVNST WYFDVPKLHP QFDYINLAAF DFLTPLRNPE EADYTAPIFF QDEQNRLPHL NVEFQVNYWL QNHCPGQKLN LGIASYGRAW KLSKESGLSG APIVHETCGA APGGIIQSA EGLLSWPEIC SKLSQNASQA YRGELAPLRK VTDLTQKYGN YALRPADENG DFGVWLSFDD PDFAGIKAVY AKSKGLGGIA LFDLSYDDFR GLCTGQKYPI VRSIKYFMG
Source	Baculovirus
Target Names	Idgf1
Protein Names	Recommended name: Chitinase-like protein Idgf1 Alternative name(s): Imaginal disk growth factor protein 1
Expression Region	21-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.