



# Recombinant Drosophila simulans Chitinase-like protein Idgf1 (Idgf1)

<b>Product Code</b>	CSB-EP839766DMJ
<b>Abbreviation</b>	Idgf1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8MX41
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila simulans (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ANNLVCYYDS TSYLRQGLAK MHTNELDLAL QFCTHLVYGY AGLKSGTLEL FSLNVDLDMF YYKDITALRQ KFPQLKILLS VGGDRDVDEA HPNKYVELLE ANRTAQQNF I DSSMILLKRN GFDGLDLAFQ LPRNKPRKVH GSLGTYWKSF KKLFTGDFVV DPLAEQHKSQ FTDLVGNIKN AFRSANLMLS LTVLPNVNST WYFDVPKLHP QFDYINLAAF DFLTPLRNPE EADYTAPIFF QDEQNRLPHL NVEFQVNYWL QNHCPGQKLN LGIASYGRAW KLSKESGLSG APIVHETCGA APGGIIQSA EGLLSWPEIC SKLSQNASQA YRGELAPLRK VTDLTQKYGN YALRPADENG DFGVWLSFDD PDFAGIKAVY AKSKGLGGIA LFDLSYDDFR GLCTGQKYPI VRSIKYFMG
<b>Source</b>	E.coli
<b>Target Names</b>	Idgf1
<b>Protein Names</b>	Recommended name: Chitinase-like protein Idgf1 Alternative name(s): Imaginal disk growth factor protein 1
<b>Expression Region</b>	21-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.