



Recombinant Human Uncharacterized protein C1orf158 (C1orf158)

Product Code	CSB-EP839778HU-B
Abbreviation	C1orf158
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N1D5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFLTAVNPQP LSTPSWQIET KYSTKVL TGN WMEERRKFTR DTDKTPQSIY RKEYIPFPDH RPDQISRWYG KRKVEGLPYK HLITHHQEPP HRYLISTYDD HYNRHGYNPG LPPLRTWNGQ KLLWLPEKSD FLLAPPTNY GLYEQLKQRQ LTPKAGLKQS TYTSSYPRPP LCAMSWREHA VPVPPHRLHP FPHF
Source	E.coli
Target Names	C1orf158
Protein Names	Recommended name: Uncharacterized protein C1orf158
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.