



# Recombinant Human WD repeat-containing protein 16 (WDR16)

<b>Product Code</b>	CSB-YP839781HU
<b>Abbreviation</b>	WDR16
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N1V2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDNKISPEAQ VAELELDAVI GFNGHVPTGL KCHPDQEHMI YPLGCTVLIQ AINTKEQNFL QGHGNNVSCL AISRSGEYIA SGQVTFMGFK ADIILWDYKN RELLARLSLH KGKIEALAFS PNDLYLVSLG GPDDGSVVVW SIAKRDAICG SPAAGLNVGN ATNVIFSRCR DEMFMTAGNG TIRVWELDLP NRKIWPTECQ TGQLKRIVMS IGVDDDDSF YLGTGGDIL KMNPRTKLLT DVGPAKDKFS LGVSAIRCLK MGGLLVGSGA GLLVFCCKSPG YKPIKIKLQ GGITSITLRG EGHQFLVGTE ESHIYRVSFT DFKETLIATC HFDAVEDIVF PFGTAELFAT CAKkdirvwh TSSNRELLRI TVPNMTCHGI DFMRDGKSII SAWNDGKIRA FAPETGRLMY VINNAHRIGV TAIATTSDCK RVISGGGEGE VRVWQIGCQT QKLEEALKEH KSSVSCIRVK RNNEECVTAS TDGTCTIIDL VRLRRNQMIL ANTLFQCVCY HPEEFQIITS GTDRKIAYWE VFDGTVIREL EGSLSGSING MDITQEGVHF VTGGNDHLVK VWDYNEGEVT HVGVGHSGNI TRIRISPGNQ YIVSVSADGA ILRWKYPYTS
<b>Source</b>	Yeast
<b>Target Names</b>	CFAP52
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 16 Alternative name(s): WD40-repeat protein up-regulated in HCC
<b>Expression Region</b>	1-620
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.