



Recombinant Human ELMO domain-containing protein 1 (ELMOD1)

Product Code	CSB-EP839788HU-B
Abbreviation	ELMOD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N336
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKHFLRMLIQ VCLYFYCKFL WRCLKFVMRK LTGRCELQRI CYNTKPGASR TMKIETSLRD SKSKLLQTSV SVHPDAIEKT IEDIMELKKI NPDVNPQLGI SLQACLLQIV GYRNLIADVE KLRREAYDSD NPQHEEMLLK LWKFLKPNT LESRISKQWC EIGFQGDDPK TDFRGMGLLG LYNLQYFAER DATAAQVLS DSLHPKCRDI TKEEISKFSK AEWEKKRMDK AIGYSFAIVG INITDLAYNL LVSGALKTHF YNIAPEAPTL SHFQQTFCYL MHEFHKFWIE EDPMDIMEFN RVREKFRKRI IKQLQNPDMA LCPHFAASEG LINM
Source	E.coli
Target Names	ELMOD1
Protein Names	Recommended name: ELMO domain-containing protein 1
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.