



Recombinant Human Uncharacterized protein C3orf23 (C3orf23)

Product Code	CSB-MP839790HU
Abbreviation	C3orf23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3R3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFCHLRPMRR LCLEKIFPHW FPFSRALSGA EAVNALRPFY FAVHPDFFGQ HPVEREINEN SLKRLSVYLE NLQKPGFKSL KPTQLTFYVR ETDQSSSDGQ EPFSTSGFRA VKFTLHTRDL LSTVLYILNS CSLSVEHIQS LNTNMHTQPL KEAKRMPDRP IKWDKSYYSF TGFKDPDEDL EQVSRVETTL TSWLDNNGKS AVKKLKNSLP LRKELDRKLD ELSHQLQLSD IRWQRSWGIA HRCSQLHSLS RLAQQNLETL KKAKGCTIIF TDRSGMSAVG HVMLGTMDVH HHWTKLFERL PSYFDLQRRRL MILEDQISYL LGGIQVVYIE ELQPVLTLEE YYSLLDVFYN RLLKSRILFH PRSLRGLQMI LNSDRYAPSL HELGHFNIPT LCDPANLQWF ILTKAQQARE NMKRKEELKV IENELIQAST KKFSLEKLYK EPSISSIQMV DCCKRLLEQS LPYLHGMHLC ISHFYSVMQD GDLCIPWNWK NGEAIK
Source	Mammalian cell
Target Names	TCAIM
Protein Names	Recommended name: Uncharacterized protein C3orf23
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.