



Recombinant Human Inactive serine protease 35 (PRSS35)

Product Code	CSB-MP839792HU
Abbreviation	PRSS35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3Z0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	LIDG SEMEWDFMWH LRKVPRIVSE RTFHILTSPAF EADAKMMVNT VCGIECQKEL PTPSLSELED YLSYETVFEN GTRTLTRVKV QDLVLEPTQN ITTKGVSVRR KRQVYGTDSR FSILDKRFLT NFPFSTAVKL STGCSGILIS PQHVLTAHC VHDGKDYVKG SKKLRVGLLK MRNKSGGKKR RGSKRSRREA SGGDQREGTR EHLRERAKGG RRRKKSGRGQ RIAEGRPSFQ WTRVKNTHIP KGWARGGMGD ATLDYDYALL ELKRAHKKKY MELGISPTIK KMPGGMIHFS GFDNDRADQL VYRFCSVSDE SNDLLYQYCD AESGSTGSGV YLRLKDPDKK NWKRKIIAVY SGHQWVDVHG VQKDYNVAVR ITPLKYAQIC LWIHGNDANC AYG
Source	Mammalian cell
Target Names	PRSS35
Protein Names	Recommended name: Inactive serine protease 35
Expression Region	17-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.