



Recombinant Human Leucine-rich repeat-containing protein 18 (LRRC18)

Product Code	CSB-BP839794HU
Abbreviation	LRRC18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N456
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVKGEKGPKG KKITLKVARN CIKITFDGKK RLDLSKMGIT TFPKCILRLS DMDLDSLRSN LIRKIPDSIS KFQNLRWLDL HSNYIDKLPE SIGQMTSLLY LNVSNRLTS NGLPVELKQL KNIRAVNLGL NHLDSVPTTL GALKELHEVG LHDNLLNNIP VSISKLPKPKL KLNIKRNPPF KPGESEIFID SIRRLNLYV VEEKDLCAAC LRKCQNARDN LNRIKNMATT TPRKTIFPNL ISPNMAKDS WEDWRIRLTS S
Source	Baculovirus
Target Names	LRRC18
Protein Names	Recommended name: Leucine-rich repeat-containing protein 18
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.