



Recombinant Human DTW domain-containing protein 1 (DTWD1)

Product Code	CSB-EP839801HU
Abbreviation	DTWD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N5C7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SLNPPIFLK RSEENSSKFV ETKQSQTTSI ASEDPLQNLC LASQEVLQKA QQSGRSKCLK CGGSRMFYCY TCYVPVENVP IEQIPLVKLP LKIDIIKHPN ETDGKSTAIH AKLLAPEFVN IYTYPCIPEY EEKDHEVALI FPGPQSISIK DISFHLQKRI QNNVRGKNDD PDKPSFKRKR TEEQEFCDLN DSKCKGTTLK KIIFIDSTWN QTNKIFTDER LQGLLQVELK TRKTCFWRHQ KGKPDFTLST IEAIYYFLVD YHTDILKEYK RGQYDNLFF YSFMYYQLIKN AKCSGDKETG KLTH
Source	E.coli
Target Names	DTWD1
Protein Names	Recommended name: DTW domain-containing protein 1
Expression Region	2-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.