



Recombinant Human Myc target protein 1 (MYCT1)

Product Code	CSB-MP839809HU
Abbreviation	MYCT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N699
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRTQVYEGLC KNYFSLAVLQ RDRIKLLFFD ILVFLSVFLL FLLFLVDIMA NNTTSLGSPW PENFWEDLIM SFTVSMAIGL VLGGFIWAVF ICLSRRRRAS APISQWSSSR RSRSSYTHGL NRTGFYRHSG CERRSNLSLA SLTFQRQASL EQANSFPRKS SFRASTFHPF LQCPPLPVET ESQLVTLPS NISPTISTSH SLSRPDYWSS NSLRVGLSTP PPPAYESI K AFPDS
Source	Mammalian cell
Target Names	MYCT1
Protein Names	Recommended name: Myc target protein 1 Alternative name(s): Myc target in myeloid cells protein 1
Expression Region	1-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.