



Recombinant Human Methyltransferase-like protein 25 (METTL25)

Product Code	CSB-MP839813HU
Abbreviation	C12orf26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N6Q8
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAASCPLPVT PDLPTLRACL QGLLQFLRDA LSISNAHTVD FYTESVWEEL VDLPPETVLA ALRKSASETE ALPSETRPLV EAEWEAGMTD FPKIFCETSQ KLVSVEAFAL AAKYYSVQNL GICTPFEQLL VALRGNQNQR IGENQKAVEF MNMKKSHEVQ AMSELISSIA DYYGIKQVID LGSGKGYLSS FLSLKYGLKV YGIDSSNTNT HGAEERNRKL KKHWWKLCHAQ SRLDVNGLAL KMAKERKVQN KVKNKADTEE VFNNSPTNQE KMPTSAILPD FSGSVISNIR NQMETHLSQP HQEENLCFEN SFSLINLLPI NAVEPTSSQQ IPNRETSEAN KERRKMTSKS SESNIYSPLT SFITADSELH DIIKDLEDCL MVGLHTCGDL APNTRLRIFTS NSEIKGVCSV GCCYHLLSEE FENQHKERTQ EKWGFPMCHY LKEERWCCGR NARMSACLAL ERVAAGQGLP TESLFYRAVL QDIIKDCYGI TKCDRHVGKI YSKCSSFLDY VRRSLKKLGL DESKLPEKII MNYEYKYPK MNELEAFNML KVLAPCIET LILLDRLCYL KEQEDIAWSA LVKLFDPVKS PRCYAVIALK KQQ
Source	Mammalian cell
Target Names	METTL25
Protein Names	Recommended name: Uncharacterized protein C12orf26
Expression Region	1-603
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.