



Recombinant Human Testis-expressed sequence 9 protein (TEX9)

Product Code	CSB-MP839814HU
Abbreviation	TEX9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N6V9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAGRSLCLTR SSVPGTPFPF PVQQPSTPGP DLLALEEEYK RLNAELQAKT ADVQQAKEI IRDRQEVRSR PVSTQMKSCD DEDDYSLRGL LPSEGIVHLH SETKPKTKNI DPNKVKQNKL HSANKGRKTN SSVKLKYSVDV QTADDVAIPE DFSDFSIAKT ISKIEGQLEE EGLPEYIDDI FSGVSNIDGT EAQIRFLKAK LHVMQEELDN VVCECNKKED EIQNLKSQVK NFEEDFMRQQ RTINMQQSQV EKYKTLFEEA NKKYDGLQQQ LSSVERELEN KRRLQKQAAS SQSATEVRLN RALEEAKEYK LELSKLRQNN KDIANEHKK IEVLKSENKK LEKQKGELMI GFKKQLKLID VLKRQKMHIE AAKMLSFTTE EFMKALEWGN S
Source	Mammalian cell
Target Names	TEX9
Protein Names	Recommended name: Testis-expressed sequence 9 protein
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.