



# Recombinant Human Putative E3 ubiquitin-protein ligase UBR7 (UBR7)

<b>Product Code</b>	CSB-BP839822HU
<b>Abbreviation</b>	UBR7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N806
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGAEGAAGR QSELEPVVSL VDVLEEDEEL ENEACAVLGG SDSEKCSYSQ GSVKRQALYA CSTCTPEGEE PAGICLACSY ECHGSHKLFY LYTKRNFRCY CGNSKFKNLE CKLLPDKAKV NSGNKYDNF FGLYCKIKRP YPDPEDEIPD EMIQCVVCEW WFHGRHLGAI PPESGDFQEM VCQACMKRCS FLWAYAAQLA VTKISTEDDG LVRNIDGIGD QEVIKPENGE HQDSTLKEDV PEQGKDDVRE VKVEQNSEPC AGSSSESDLQ TVFKNESLNA ESKSGCKLQE LKAKQLIKKD TATYWPLNWR SKLCTCQDCM KMYGDLVDLF LTDEYDVLVLA YENKGKIAQA TDRSDPLMDT LSSMNRVQQV ELICEYNDLK TELKDYLKRF ADEGTVVKRE DIQQFFEEFQ SKKRRRVDGM QYYCS
<b>Source</b>	Baculovirus
<b>Target Names</b>	UBR7
<b>Protein Names</b>	Recommended name: Putative E3 ubiquitin-protein ligase UBR7 EC= 6.3.2.- Alternative name(s): N-recognin-7
<b>Expression Region</b>	1-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.