



Recombinant Human Uncharacterized protein C7orf31 (C7orf31)

Product Code	CSB-MP839824HU
Abbreviation	C7orf31
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N865
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVIHGRPYC CRELEGADIL SNTFYSNELH NPLQTVTRPT ASEDYQELR ESLQQCRLPW GAEREYGGII PISLPEDHRP KYEPPRVMGK GHQHYGFGGE TWPRKLPVEQ FYYLTQNKKS DVYGNDSLIP KPPNSTVGEI CLPYPIEHPY HTHICRGAMF PTFTSPEDLY TGIKARTQQP FPPTVPTKAY DSTVLKTRGN PYRYELIDIP MDSKKAALTW PGQGVYYDFP RGVEKNKPVF YPKPPKTFAP NTSLNSWDPI CSAKEANIQR NLERSHWLTS YTHDFTGLGP MDPLELDDYH EKMVAELTRK IGFDPEPQEK FHPVFKPPRP LEGRIARLIQ NRRSLEAIVQ QRPRSCPDCT PRVLCNFHTF VPSSKEMVAL SDNIPAGVTH KNQDIEEKII EEQSLLSTYE LPSCYPTKDL TSIYDIKPPF KITDTKKTED LYWRQQSLKT QPTPYCKPDH WIHYENLKSP LRDQYNMCPD PVSLSKPSVL QNKQDTEAFT LEHFLSKPEE ELFLNMENNE ETRPVLGWIP RAGVTKPQTN LLELKNSFSK TGAQKRFHKS ILEDHKDLRD NEHSGMKHQF YGHNSYFFYN
Source	Mammalian cell
Target Names	C7orf31
Protein Names	Recommended name: Uncharacterized protein C7orf31
Expression Region	1-590
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.