



Recombinant Human Protein deltex-3 (DTX3)

Product Code	CSB-MP839833HU
Abbreviation	DTX3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8N9I9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSFVLSRMAA CGGTCKNKVT VSKPVWDFLS KETPARLARL REEHRVSILI DGETSDIYVL QLSPQGPPPA PPNGLYLARK ALKGLLKEAE KELKKAQRQG ELMGCLALGG GGEHPEMHRA GPPPLRAAPL LPPGARGLPP PPPPLPPPLP PRLREEAEEQ ESTCPICLGE IQNAKTLEKC RHSFCEGCIT RALQVKKACP MCGRFYGQLV GNQPQNGRML VSKDATLLLP SYEKYGTIVI QYVFPPGVQG AEHPNPGVRY PGTTRVAYLP DCPEGNKVLT LFRKAFDQRL TFTIGTSMTT GRPNVITWND IHHKTSCTGG PQLFGYPDPT YLTRVQEELR AKGITDD
Source	Mammalian cell
Target Names	DTX3
Protein Names	Recommended name: Protein deltex-3 Short name= Deltex3 Alternative name(s): RING finger protein 154
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>