



Recombinant Human [Protein ADP-ribosylarginine] hydrolase-like protein 1 (ADPRHL1)

Product Code	CSB-EP839860HU-B
Abbreviation	ADPRHL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NDY3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEKFKAAMLL GSVGDALGYR NVCKENSTVG MKIQEELQRS GGLDHLVLSP GEWPVSDNTI MHIATAEALT TDYWCLDDLY REMVRCYVEI VEKLPERRPD PATIEGCAQL KPNNYLLAWH TPFNEKGSGF GAATKAMCIG LRYWKPERLE TLIEVSVECG RMTHNHPTGF LGSLECTALFV SFAAQGKPLV QWGRDMLRAV PLAEEYCRKT IRHTAEYQEH WFYFEAKWQF YLEERKISKD SENKAIFPDN YDAEEREKTY RKWSSEGRGG RRGHDAPMIA YDALLAAGNS WTELCHRAMF HGGESAATGT IAGCLFGLLY GLDLVPKGLY QDLEDKEKLE DLGAALYRLS TEEK
Source	E.coli
Target Names	ADPRHL1
Protein Names	Recommended name: [Protein ADP-ribosylarginine] hydrolase-like protein 1 Alternative name(s): ADP-ribosylhydrolase 2
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.