



Recombinant Human Serum response factor-binding protein 1 (SRFBP1)

Product Code	CSB-YP839864HU
Abbreviation	SRFBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEF9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AQPGTLNLN NEVVKMRKEV KRIRVLVIRK LVRSVGRLKS KKGTEDEALLK NQRRAQRLL EIHAMKELKP DIVTKSALGD DINFEKIFKK PDSTATERAI ARLAVHPLLK KKIDVLKAAV QAFKEARQNV AEVESSKNAS EDNHSENTLY SNDNGSNLQR EATVISEQKV KETKILAKKP IHNSKEKIAK MEHGPKAVTI ANSPSKPSEK DSVVSLESQK TPADPKLCTL SQTCKNKGSD SSLSGNSDGG EEFCEEEKEY FDDSTEERFY KQSSMSSEDSD SGDDFFIGKV RRTRKKESSC HSSVKEQKPL EKVFLKEDTG ETHGDTRNDK IKPSTETRKL ESVFFHSLSG SKSSRRNFKE QAPKTRSLDF PQNEPQIKNQ FNKKLSGRLE NTKQQLQLPL HPSWEASRRR KEQQSNIAVF QGKKITFDD
Source	Yeast
Target Names	SRFBP1
Protein Names	Recommended name: Serum response factor-binding protein 1 Alternative name(s): SRF-dependent transcription regulation-associated protein p49/STRAP
Expression Region	2-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.