



Recombinant Human Beta-1,3-N-acetylglucosaminyltransferase lunatic fringe (LFNG), Yeast, partial, Yeast

Product Code	CSB-YP839868HU
Abbreviation	LFNG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NES3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	LFNG
Protein Names	Recommended name: Beta-1,3-N-acetylglucosaminyltransferase lunatic fringe EC= 2.4.1.222 Alternative name(s): O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	full length of Isoform 2
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.