



Recombinant Human Uncharacterized protein C15orf43 (C15orf43)

Product Code	CSB-EP839917HU-B
Abbreviation	C15orf43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NHR7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFQGQRGWFC GSVSQDLRQF WVAEGGTISD PRAADFLFSC DASHPDTLRI YQSLDYIEDN ATVFHAYYLS AVANAKIKNS VALGHFILPP ACLQKEIRRK IGSFIWEQDQ HFLIEKHDEV TPNEIKTLRE NSELATEHKK ELSKSPEKHF IRTPVVEKQM YFPLQNYPVN NMVTGYISID AMKKFLGELH DFIPGTSGYL AYHVQNEINM SAIKNKLKRK
Source	E.coli
Target Names	TERB2
Protein Names	Recommended name: Uncharacterized protein C15orf43
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.