



Recombinant Arabidopsis thaliana Cyclin-L1-1 (CYCL1-1)

Product Code	CSB-MP840426DOA
Abbreviation	CYCL1-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8RWW3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MIYTAIDNFY LSDEQLKASP SRKDGIDETT EISLRIYGCD LIQEGGILLK LPQAVMATGQ VLFQRFYCKK SLAKFDVKIV AASCVWLASK LEENPKKARQ VIIVFHRMEC RRENLPLEHL DMYAKKFSEL KVLSRTERH ILKEMGFVCH VEHPHKFISN YLATLETPPE LRQEAWNLAN DSLRRTLCVR FRSEVVACGV VYAAARRFQV PLPENPPWWK AFDADKSSID EVCRLAHL SLPKAQYISV CKDGKPFDFS SRSGNSQGQS ATKDLLPGAG EAVDTKCTAG SANNDLKDGM VTTPEKATD SKKSGTESNS QPIVGDSSYE RSKVGDRERE SDREKERGRE RDRGRSHRGR DSDRDSDRER DKLKDRSHHR SRDRLKDSGG HSDKSRHHSS RDRDYRDSSK DRRRHH
Source	Mammalian cell
Target Names	CYCL1-1
Protein Names	Recommended name: Cyclin-L1-1 Short name= Cycl1;1
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.