



Recombinant *Oryza sativa* subsp. japonica Thioredoxin F, chloroplastic (Os01g0913000, LOC_Os01g68480)

Product Code	CSB-MP840439OFG
Abbreviation	Os01g0913000, LOC_Os01g68480
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8S091
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AEAVAVTG QVTEVNKDTF WPIVKSAGPK VVVLDMYTQW CGPCKVMAPK FQEMSEKDQD VVFLKLDCNQ DNKSLAKELG IKVVPTFKIL KDGKVVKEVT GAKLDELIQA IETVKSS
Source	Mammalian cell
Target Names	Os01g0913000
Protein Names	Recommended name: Thioredoxin F, chloroplastic Short name= OsTrxf Alternative name(s): OsTrx03
Expression Region	73-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.