



Recombinant Arabidopsis thaliana Solanesyl diphosphate synthase 1 (SPS1)

| | |
|--------------------------|--|
| Product Code | CSB-EP840457DOA-B |
| Abbreviation | SPS1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8S948 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MMTSCRNIDL GTMMMACGCG RRQFPSLAKT VCKFTSSNRS YGGLVGSCKA VPTKSKEISL LNGIGQSQTV SFDLKQESKQ PISLVTLFEL VAVDLQTLND NLLSIVGAEN PVLISAAEQI FGAGGKMRMP GLVFLVSHAT AELAGLKELT TEHRLAEII EMIHTASLIH DDVLDESDMR RGKETVHELFF GTRVAVLAGD FMFAQASWYL ANLENLEVIK LISQVIK DFA SGEIKQASSL FDCDTKLDEY LLKSFYKTAS LVAASTKGAA IFSRVEPDVT EQMYEFGKNL GLSFQIVDDI LDFTQSTEQL GKPAGSDLAK GNLTAPVIFA LEREPRLEI IESEFCEAGS LEEAIEAVTK GGGIKRAQEL AREKADDAIK NLQCLPRSGF RSALEDMVLY NLERID |
| Source | E.coli |
| Target Names | SPS1 |
| Protein Names | Recommended name: Solanesyl diphosphate synthase 1 Short name= AtSPS1 EC= 2.5.1.85 Alternative name(s): All-trans-nonaprenyl-diphosphate synthase 1 (geranylgeranyl-diphosphate specific) |
| Expression Region | 1-406 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.