



Recombinant *Oryza sativa* subsp. japonica 50S ribosomal protein L18, chloroplastic (RPL18)

Product Code	CSB-YP840460OFG
Abbreviation	RPL18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SAY0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	LRKKVSG TTERPRLSVF RSNKHLYAQV IDDTKSCTLV SASTMHKSLK KDLEYSAGPT VEVAQKIGEV IAKSCLEKGI TKVVFDGGF LYHGRIKALA DAARENLDF
Source	Yeast
Target Names	RPL18
Protein Names	Recommended name: 50S ribosomal protein L18, chloroplastic Alternative name(s): CL18
Expression Region	64-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.