



Recombinant Encephalitozoon cuniculi Probable cell division control protein 7 homolog 1 (CDC7-1)

Product Code	CSB-BP840489EKH
Abbreviation	CDC7-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SR85
Product Type	Recombinant Protein
Immunogen Species	Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)
Purity	>85% (SDS-PAGE)
Sequence	MEEKILESD LRHISFVMPK YTPIEKIGEG SFSVVYKALD AESGRYVALK AITRTSSPAR VLDEMMFLKT LGGRKNCMGL LGCFRNEDQV VAVFPYFPEI DFREFISNAN LADIKRYLHN LLIAIEHVHS NGIMHRDLKP GNFLYNKESG RGMLIDFGLA QYEEYSEGQH AEGGAKPAGP LFFNSVVSK TKPPGYERD GRPPMKAPRA GTRGFRAPEV LFRCQRQTGA IDMWSVGVIF LTILTTQYPF FYSSDDIDS IVEIATIFGHA EMRKAAKFYG RVWRSNIDS IPEERIPFETI VESLNPWAEI GSDGYDLLYR MLDLCSSSRI T ASDALSHPF FDDLKTHENC A
Source	Baculovirus
Target Names	CDC7-1
Protein Names	Recommended name: Probable cell division control protein 7 homolog 1 EC= 2.7.11.1
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.