



Recombinant Encephalitozoon cuniculi Probable serine/threonine-protein kinase KIN28 homolog (KIN28)

Product Code	CSB-MP840527EKH
Abbreviation	KIN28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SW92
Product Type	Recombinant Protein
Immunogen Species	Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)
Purity	>85% (SDS-PAGE)
Sequence	MKTYIRERRL GEGTYAVIYL GYRALPQDKP LVSSGTRIED VPVAIKKIKP TKYTQGHEIS AIREIKSLKR IDSKYVVRLL DTFVYDKCVH IVLEYVETNL ENVIRNSDKI IMPGDIKAWI LMVLRGVYEC HRLFIIHRDI KPNNILITSE GMVKLADFGL TRGIGNRMTP QAVTRWYRAP ELLMGSRDYG SPVDMWSVGC VFAELFLRVP LFAGDTDIQQ LDMIFRALGT PVEREWPGVS TLPEFLDFQQ YPKASLNGLF SAASSDALDL LEKLLTLNPC NRISCDDAIK HPYFKSSPPP TPIGKLPV
Source	Mammalian cell
Target Names	KIN28
Protein Names	Recommended name: Probable serine/threonine-protein kinase KIN28 homolog EC= 2.7.11.23
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.