



Recombinant Human Protein FAM81A (FAM81A)

Product Code	CSB-MP840559HU
Abbreviation	FAM81A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TBF8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MENMHLRRVR TMPRHSQSLT MAPYSSVSLV EQLEDRLCH EKTTAALVEH AFRIKDDIVN SLQKMQNKG GDRLARLFLE EHIRNITAIV KQLNRDIEVL QEQRARDNI SYGTNSALKT LEMRQLSGLG DLRGRVARCD ASIARLSAEH KTTYEGLQHL NKEQQA AKLI LETKIKDAEG QISQLLNVRD LSISEQSTKL KMSHRDSNHQ LQLLDTKFKG TVEELSNQIL SARSWLQQEQ ERIEKELLQK IDQLSLIVKE NSGASERDME KKLSQMSARL DKIEEGQKKT FDGQRTRQEE EKMHGRTKL ELQMNQNIKE MKAENVAGFT AVYESIGSLR QVLEAKMKLD RDQLQKQIQL MQKPETPM
Source	Mammalian cell
Target Names	FAM81A
Protein Names	Recommended name: Protein FAM81A
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.