



Recombinant Human Protein FAM45A (FAM45A)

Product Code	CSB-BP840565HU
Abbreviation	FAM45A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TCE6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAEVADTQ LMLGVGLIEK DTNGEVLWVW CYPSTTATLR NLLLRKCCLT DENKLLHPFV FGQYRRTWFY ITTIEVPDSS ILKKVTHFSI VLTAKDFNPE KYAAFTRILC RMYLKHGSPV KMMESYIAVL TKGICQSEEN GSFLSKDFDA RKAYLAGSIK DIVSQFGMET VILHTALMLK KRIVVYHPKI EAVQEFTRTL PALVWHRQDW TILHSYVHLN ADELEALQMC TGYVAGFVDL EVSNRPDLVD VFNLAESEI TIAPLAKEAM AMGKLHKEMG QLIVQSAEDP EKSESHVIQD IALKTREIFT NLAPFSEVSA DGEKRVLNLE ALKQKRFFPA TENFLYHLAA AEQMLKI
Source	Baculovirus
Target Names	FAM45A
Protein Names	Recommended name: Protein FAM45A
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.