



Recombinant Human BPI fold-containing family B member 1 (BPIFB1)

Product Code	CSB-BP840572HU
Abbreviation	BPIFB1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDL5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> TLSPTAVLI LGPKVIKEKL TQELKDHNAT SILQQLPLLS AMREKPAGGI PVLGSLVNTV LKHIIWLKVI TANILQLQVK PSANDQELLV KIPLDMVAGF NTPLVKTIVE FHMTTEAQAT IRMDTSASGP TRLVLSDCAT SHGSLRIQLL HKLSFLVNAL AKQVMNLLVP SLPNLVKNQL CPVIEASFNG MYADLLQLVK VPISLSIDRL EFDLLYPAIK GDTIQLYLGA KLLDSQGKVT KWFNNSAASL TMPTLDNIPF SLIVSQDVVK AAVAAVLSPE EFMVLLDSVL PESAHRLKSS IGLINEKAAD KLGSTQIVKI LTQDTPEFFI DQGHAKVAQL IVLEVFPSSE ALRPLFTLGI EASSEAQFYT KGDQLILNLN NISSDRIQLM NSGIGWFQPD VLKNIITEII HSILLPNQNG KLRSGVPVSL VKALGFEEAE SSLTKDALVL TPASLWKPSS PVSQ </p>
Source	Baculovirus
Target Names	BPIFB1
Protein Names	Recommended name: BPI fold-containing family B member 1 Alternative name(s): Long palate, lung and nasal epithelium carcinoma-associated protein 1 Von Ebner minor salivary gland protein Short name= VEMSGP
Expression Region	22-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.