



Recombinant Human U3 small nucleolar RNA-associated protein 15 homolog (UTP15)

Product Code	CSB-BP840580HU
Abbreviation	UTP15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TED0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AGYKPVAIQ TYPILGEKIT QDTLYWNNYK TPVQIKEFGA VSKVDFSPQP PYNAYVTASS RIHIYGRYSQ EPIKTFSRFK DTAYCATFRQ DGRLLVAGSE DGGVQLFDIS GRAPLRQFEG HTKAVHTVDF TADKYHVVSG ADDYTVKLWD IPNSKEILTF KEHSDYVRCG CASKLNPDFL ITGSYDHTVK MFDARTSESV LSVEHGQPVE SVLLFPSGGL LVSAGGRYVK VWDMLKGGQL LVSLKNHHKT VTCLCLSSSG QRLLSGSLDR KVKVYSTTSY KVVHSFDYAA SILSLALAHE DETIVVGMTN GILSVKHRKS EAKKESLPRR RRPAYRTFIK GKNYMKQRDD ILINRPAKKH LELYDRDLKH FRISKALDRV LDPTCTIKTP EITVSIKEL NRRGVLANAL AGRDEKEISH VLNFLIRNLS QPRFAPVLIN AAEIIDIYL PVIGQSPVVD KKFLLLQGLV EKEIDYQREL LETLGMMMDML FATMRRKEGT SVLEHTSDGF PENKKIES
Source	Baculovirus
Target Names	UTP15
Protein Names	Recommended name: U3 small nucleolar RNA-associated protein 15 homolog
Expression Region	2-518
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.