



Recombinant *Pyrococcus furiosus* DNA double-strand break repair protein Mre11 (mre11)

Product Code	CSB-BP840728FHW
Abbreviation	mre11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8U1N9
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrococcus furiosus</i> (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1)
Purity	>85% (SDS-PAGE)
Sequence	MKFAHLADIH LGYEQFHKPQ REEF AEAFK NALEIAVQEN VDFILIAGDL FHSSRPSPGT LKKAIALLLQI PKEHSIPVFA IEGNHDRTRQ GPSVLNLLD FGLVYVIGMR KEKVENEYLT SERLNGEYL VKGVYKDLEI HGMKYMSSAW FEANKEILKR LFRPTDNAIL MLHQGVREVS EARGEDYFEI GLGDLPEGYL YYALGHIHKR YETSYSGSPV VYPGSLERWD FGDYEVRYEW DGIKFKERYG VNKGFYIVED FKPRFVEIKV RPFIDVKIKG SEEEIRKAIK RLIPLIPKNA YVRLNIGWRK PFDLTEIKEL LNVEYLKIDT WRIKERTDEE SGKIGLPSDF FTEFELKIID ILGEKDFDDF DYIIKLITEG KVEEEGPL EE AVKKVSEEKG KTVRQKIESI PKKKRGTLD S WLGGAR
Source	Baculovirus
Target Names	mre11
Protein Names	Recommended name: DNA double-strand break repair protein Mre11 Alternative name(s): pfMre11
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.