



Recombinant Mouse Aldo-keto reductase family 1 member C13 (Akr1c13)

Product Code	CSB-MP840826MO
Abbreviation	Akr1c13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VC28
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSSKQHCVKL NDGHLIPALG FGTYKPKEVP KSKSLEAACL ALDVGYRHVD TAYAYQVEEE IGQAIQSKIK AGVVKREDLF ITTKLWCTCF RPELVKPALE KSLKKLQLDY VDLYIMHYPV PMKSGDNDFP VNEQGKSLLD TVDFCDTWER LEECKDAGLV KSIGVSNFNH RQLERILNKP GLKYKPCVNQ VECHLYLNQR KLLDYCESKD IVLVAYGALG TQRYKEWVDQ NSPVLNNDPV LCDVAKKNKR SPALIALRYL IQRGIVPLAQ SFKENEMREN LQVFGFQLSP EDMKTLDGLN KNFRYLPAEF LVDHPEYPFV EEY
Source	Mammalian cell
Target Names	Akr1c13
Protein Names	Recommended name: Aldo-keto reductase family 1 member C13 EC= 1.1.1.-
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.