



Recombinant Mouse E3 ubiquitin-protein ligase RNF8 (Rnf8)

Product Code	CSB-BP840827MO
Abbreviation	Rnf8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VC56
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGEPDPLVSG QLAARRSWCL RRLGMDCEWL QLEAGTEVTI GRGLSVTYQL ISKVCPLMIS RSHCVLKQNP EGQWTIMDNK SLNGVWLNRE RLAPLQGYCI RKGDHIQLGV PLESRETAEY EYEVIEEDWE SLAPCLAPKN DQRMEKHKGS RTKRKFSSPG LENLPAEGSS DLRCPLANVA SKPIEPEKLH GKGDASSQSL GCLCPGLTSL KASERAAGPH ACSALPKVLE LSCPCKKQKAC RPSASQNSLE LFKVTMSRML KLKTQMKEKQ IAVLNVKRQT RKGSSKKIVR MEKELRNLQS QLYAEQAQQQ ARVEQLEKTF QEEAHYLQGL EKEQGECDLK QQLVQALQEH QALMEELNCS KKDFEKIIQA KNKELEQTKE EKDKVQAQKE EVLSHMNDLL ENELQCIICS EYFIEAVTLN CAHSFCSFCI NEWMKRKVEC PICRKDIESR TNSLVLDNCI SKMVDNLSSD VKERRSVLIR ERRAKRLS
Source	Baculovirus
Target Names	Rnf8
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF8 EC= 6.3.2.- Alternative name(s): ActA-interacting protein 37 Short name= AIP37 LaXp180 RING finger protein 8
Expression Region	1-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein contains a RING finger motif and a FHA domain. This protein has been shown to interact with several class II ubiquitin-conjugating enzymes (E2), including UBE2E1/UBCH6, UBE2E2, and UBE2E3, and may act as an ubiquitin ligase (E3) in the ubiquitination of certain nuclear proteins. Alternatively spliced transcript variants encoding distinct isoforms have been reported.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.