



Recombinant Bovine Protein S100-A9 (S100A9)

| | |
|--------------------------|---|
| Product Code | CSB-EP020642BO-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P28783 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | EDKMSQMES SIETIINIFH QYSVRLGHYD TLIQKEFKQL VQKELPNFLK KQKKNEAAIN EIMEDLDTNV DKQLSFEEFI MLVARLTVAS HEEMHNTAPP GQGHRHGPGY GKGGSGSCSG QGSPDQGSHD LGSHGHGHGH SHGGHGHSHG GHGSH |
| Source | E.coli |
| Target Names | S100A9 |
| Protein Names | Recommended name: Protein S100-A9 Alternative name(s): BEE22 Calgranulin-B Neutrophil cytosolic 23 kDa protein Short name= p23 S100 calcium-binding protein A9 |
| Expression Region | 2-156 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This protein is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in the inhibition of casein kinase and altered expression of this protein is associated with the disease cystic fibrosis. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |