



Recombinant Rat FERM domain-containing protein 6 (Frmd6)

Product Code	CSB-MP840879RA
Abbreviation	Frmd6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VII0
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MNKLTFHNNK VMQDRRRVCI FLPNDKSVSI IINVKILCHQ LLVQVCDLLR LKDSLHFLGLS VIQNNEHVYM ELSQKLYKYC PKEWKKEASK GIDQFGPPMI IHFRVQYYVE NGKLISDRIA RYYYYWHLRK QVLHSQCCLR EEAYFLLAAF ALQADLGNFK RKVHHGDYFE PEAYFPAWVV SKRGKDYILK HIPNMHRDQF ALTASEAYLK YIKEAVRLDD VAIHYYRLYK DKREVEGSLT LGLTMRGIQI FQNLEEEKQL LYDFPWTNNG KLVFVGKKFE ILPDGLPSAR KLVYYTGCPT RSRHLQLLS NSHRLYMNLQ PVLRLHLR
Source	Mammalian cell
Target Names	Frmd6
Protein Names	Recommended name: FERM domain-containing protein 6 Alternative name(s): FERM domain-containing protein 163SCII
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.