



Recombinant Mouse Nitrogen permease regulator 3-like protein (Nprl3)

Product Code	CSB-BP840880MO
Abbreviation	Nprl3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VIJ8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MGDNTSPISV ILVSSGSRGN KLLFRYPFQR SQEHPASQTN KPRSRVAVNN TGEHADDQDG DSRFSDVILA TILATKSEMC GQKFELKIDN VRFVGHPTLL QHALGQVSKT DPSPKREAPT MILFNVVFAL RANADPSVIN CLHNLSRRIA TVLQHEERRC QYL TREAKLI LALQDEVSAM ADANEGPQSP FQHILPKCKL ARDLKEAYDS LCTSGVVR LH INSWLEVSFC LPHKIHYAAS SLIPPEAIER SLKAIRPYHA LLLL SDEKSL LSELPIDCSP ALVRVIKTT S AVKNLQQLAQ DADLALLQVF QLA AHLVYWG KAVIYPLCE NNVYVMSPNA SVCLYSPLAE QFSRQFPSHD LPSVLAKFSL PVSLSEFRSP LAPPAQETQL IQMVVWMLQR RLLIQLHTYV CLMASPSEEE PRLREDDVPF TARVGGRSLS TPNALSFGSP TSSDDMTLTS PSMDNSSAEL LPSGDSPLNK RMTENLLASL SEHERAAILN VPAAQNPEDL RMFARLLHYF RGRHHLEEIM YNENTRRSQL LMLFDKFRSV LVVTTHEDPV IAVFQALLT
Source	Baculovirus
Target Names	Nprl3
Protein Names	Recommended name: Nitrogen permease regulator 3-like protein Alternative name(s): Alpha-globin regulatory element-containing gene protein
Expression Region	1-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.