



Recombinant Arabidopsis thaliana Pyridoxal kinase (PK)

Product Code	CSB-BP840922DOA
Abbreviation	PK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W1X2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTTPPVLSLA LPSDTGRVLS IQSHTVQGYV GNKSAVFPLQ LLGYDVDPIN SVQFSNHTGY PTFKGQVLNG QQLCDLIEGL EANDLLFYTH VLTGYIGSVS FLDTILEVIN KLRSVNPNTL YVCDPVMGDE GKLYVPEELV HVYREKVVPL ASMLTPNQFE AEKLTGLRIN SEEDGREACA ILHAAGPSKV VITSITIGGI LLIGSHQKE KGLKPEQFKI LIHKIPAYFT GTGDLMTALL LGWSNKYPD LDKAAELAVS TLQALLRRTL DDYKRAGYDP TSSSLEIRLI QSQEDIRNPK VELKAERYS
Source	Baculovirus
Target Names	PK
Protein Names	Recommended name: Pyridoxal kinase EC= 2.7.1.35 Alternative name(s): Protein SALT OVERLY SENSITIVE 4 Pyridoxal kinase-like protein SOS4 Pyridoxine kinase
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.