



Recombinant Arabidopsis thaliana Putative lipid-transfer protein DIR1 (DIR1)

Product Code	CSB-MP840926DOA
Abbreviation	DIR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W453
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	AIDLC GMSQDELNEC KPAVSKENPT SPSQPCCTAL QHADFACLGG YKNPWLGSF GVDPELASAL PKQCGLANAP TC
Source	Mammalian cell
Target Names	DIR1
Protein Names	Recommended name: Putative lipid-transfer protein DIR1 Alternative name(s): Protein DEFECTIVE IN INDUCED RESISTANCE 1
Expression Region	26-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.