



# Recombinant Human Sec1 family domain-containing protein 2 (SCFD2)

<b>Product Code</b>	CSB-EP840978HU
<b>Abbreviation</b>	SCFD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WU76
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSASGVLSTF QQGWEQVLAK VKRAVVYLDA ACAESLHWGC GSTRLLLEAVG GPDCHLREFE PDAIGGGAKQ PKAVFVLSCL LKGRTVEILR DIICRSHFQY CVVVTTVSHA VHILTANHVPA AAAAEMEGQQ PVFEQLEEKL CEWMGNMNYT AEFVHVPLLL APVAPHFALT PAFASLFPLL PQDVHLLNSA RPDKRKLGS L GDVDSTTLTP ELLLQIRCLV SGLSSLCEHL GVREECFAVG SLSQVIAADL ANYAPAKNRK KTAAGRASVV FVDRTLDTLG AVGHHGDNLV EKIISALPQL PGHTNDVMVN MIALTALHTE EENYNVVAPG CLSQSSDTTA KALWEALLNT KHKEAVMEVR RHLVEAASRE NLPIKMSMGR VTPGQLMSYI QLFKNNLKAL MNHCGLLQLG LATAQTLKHP QTAKWDNFLA FERLLLQSIG ESAMSVVLNQ LLPMIKPVTQ RTNEDYSPEE LLLIIYYS VTGELTVDKD LCEAEEKVKK ALAQVFCEES GLSPLLQKIT DWDSSINLTF HKSKIAVDEL FTSLRDIAGA RSLLKQFKSV YVPGNHTHQA SYKPLLKQVV EEIFHPPERPD SVDIEHMSSG LTDLLKTGFS MFMKVS RPHP SDYPLLLIFV VGGVTVSEVK MVKDLVASLK PGTQVIVLST RLLKPLNIPE LLFATDRLHP DLGF
<b>Source</b>	E.coli
<b>Target Names</b>	SCFD2
<b>Protein Names</b>	Recommended name: Sec1 family domain-containing protein 2 Alternative name(s): Syntaxin-binding protein 1-like 1
<b>Expression Region</b>	1-684
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.