



Recombinant Human Sec1 family domain-containing protein 2 (SCFD2)

Product Code	CSB-MP840978HU
Abbreviation	SCFD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WU76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSASGVLSTF QQGWEQVLAK VKRAVVYLDA ACAESLHWGC GSTRLLLEAVG GPDCHLREFE PDAIGGGAKQ PKAVFVLSCL LKGRTVEILR DIICRSHFQY CVVVTTVSHA VHILTANHVPA AAAAEMEGQQ PVFEQLEEKL CEWMGMNMYT AEFVHVPLLL APVAPHFALT PAFASLFPLL PQDVHLLNSA RPDKRKLGS L GDVDSTTLTP ELLLQIRCLV SGLSSLCEHL GVREECFAVG SLSQVIAADL ANYAPAKNRK KTAAGRASVV FVDRTLDTLG AVGHHGDNLV EKIISALPQL PGHTNDVMVN MIALTALHTE EENYNVVAPG CLSQSSDTTA KALWEALLNT KHKEAVMEVR RHLVEAASRE NLPIKMSMGR VTPGQLMSYI QLFKNNLKAL MNHCGLLQLG LATAQTLKHP QTAKWDNFLA FERLLLQSIG ESAMSVVLNQ LLPMIKPVTQ RTNEDYSPEE LLLIYIYS VTGELTVDKD LCEAEEKVKK ALAQVFCEES GLSPLLQKIT DWDSSINLTF HKSKIAVDEL FTSLRDIAGA RSLKQFKSV YVPGNHTHQA SYKPLLKQVV EEIFHPPERPD SVDIEHMSSG LTDLLKTGFS MFMKVS RPHP SDYPLLLIFV VGGVTVSEVK MVKDLVASLK PGTQVIVLST RLLKPLNIPE LLFATDRLHP DLGF
Source	Mammalian cell
Target Names	SCFD2
Protein Names	Recommended name: Sec1 family domain-containing protein 2 Alternative name(s): Syntaxin-binding protein 1-like 1
Expression Region	1-684
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.