



Recombinant Human Uncharacterized protein C19orf29 (C19orf29)

Product Code	CSB-BP840982HU
Abbreviation	C19orf29
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUQ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGRDTRSRSR SAGRRGRRRQ SQSGSRSRSR SHGRRNRRRR EDEGRRRRRR RSRERRSDSE EERWQRSGMR SRSPRPKWH SRDGSSQSDS GEEQSRGQWA RRRRRARSWS PSSSASSAS PGRSQSPRAA AAALSQQQSL QERLRLREER KQQEELMKAF ETPEEKRRARR LAKKEAKERK KREKMGWGEE YMGYTNTDNP FGDNNLLGTF IWNKALEKKG ISHLEEKELK ERNKRIQEDN RLELQKVKQL RLEREREKAM REQELEMLQR EKEAEHFKTW EEQEDNFHLQ QAKLRSKIRI RDGRAKPIDL LAKYISAEDD DLAVEMHEPY TFLNGLTVAD MEDLLEDIQV YMELEQGKNA DFWRDMMTTIT EDEISKLRKL EASGKGPGER REGVNASVSS DVQSVFKGKT YNQLQVIFQG IEGKIRAGGP NLDMGYWESL LQQLRAHMAR ARLRERHQDV LRQKLYKLLKQ EQGVESEPLF PILKQEPQSP SRSLEPEDAA PTPPGPSSEG GPAAEAVDGA TPTEGDGDGD GEGEGEGEAV LMEEDLIQQS LDDYDAGRYS PRLTAHELP LDAHVLEPDE DLQRLQLSRQ QLQVTGDASE SAEDIFFRRA KEGMGQDEAQ FSVEMPLTGK AYLWADKYRP RKPRFFNRVH TGFEWNKYNQ THYDFDNPPP KIVQGYKFNI FYPDLIDKRS TPEYFLEACA DNKDFAILRF HAGPPYEDIA FKIVNREWEY SHRHGFRQCQF ANGIFQLWFH FKRYRYRR
Source	Baculovirus
Target Names	CACTIN
Protein Names	Recommended name: Uncharacterized protein C19orf29 Alternative name(s): Cactin Renal carcinoma antigen NY-REN-24
Expression Region	1-758
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.