



# Recombinant Human Uncharacterized protein C19orf29 (C19orf29)

<b>Product Code</b>	CSB-MP840982HU
<b>Abbreviation</b>	C19orf29
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WUQ7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGRDTRSRSR SAGRRGRRRQ SQSGSRSRSR SHGRRNRRRR EDEGRRRRRR RSRERRSDSE EERWQRSGMR SRSPRPKWH SRDSSQSDS GEEQSRGQWA RRRRRARSWS PSSSASSAS PGRSQSPRAA AAALSQQQSL QERLRLREER KQQEELMKAF ETPEEKRRR LAKKEAKERK KREKMGWGEE YMGYTNTDNP FGDNNLLGTF IWNKALEKKG ISHLEEKELK ERNKRIQEDN RLELQKVKQL RLEREREKAM REQELEMLQR EKEAEHFKTW EEQEDNFHLQ QAKLRSKIRI RDGRAKPIDL LAKYISAEDD DLAVEMHEPY TFLNGLTVAD MEDLLEDIQV YMELEQGKNA DFWRDMMTTIT EDEISLRLKL EASGKGPGER REGVNASVSS DVQSVFKGKT YNQLQVIFQG IEGKIRAGGP NLDMGYWESL LQQLRAHMAR ARLRERHQDV LRQKLYKLLQ EQGVESEPLF PILKQEPQSP SRSLEPEDAA PTPPGPSSEG GPAEAEVDGA TPTEGDGDGD GEGEGEGEAV LMEEDLIQQS LDDYDAGRYS PRLTAHELP LDAHVLEPDE DLQRLQLSRQ QLQVTGDASE SAEDIFFRRA KEGMGQDEAQ FSVEMPLTGK AYLWADKYRP RKPRFFNRVH TGFEWNKYNQ THYDFDNPPP KIVQGYKFNI FYPDLIDKRS TPEYFLEACA DNKDFAILRF HAGPPYEDIA FKIVNREWEY SHRHGFRQCQF ANGIFQLWFH FKRYRYRR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CACTIN
<b>Protein Names</b>	Recommended name: Uncharacterized protein C19orf29 Alternative name(s): Cactin Renal carcinoma antigen NY-REN-24
<b>Expression Region</b>	1-758
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.