



Recombinant Human Cytoskeleton-associated protein 2 (CKAP2)

Product Code	CSB-MP840993HU
Abbreviation	CKAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWK9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSTPAVPQDL QLPPSQRAQS AFKEQRRQKL KEHLLRRKTL FAYKQENEML SSSRDQRVVT SEDQVQEGTK VLKCLKTKMAD KENMKRPAES KNNTVVGKHC IPLKPSNELT NSTVVIDTHK PKDSNQTPHL LLTEDDPQSQ HMTLSQAFHL KNNSSKKKQMT TEKQKQDANM PKKPVLGSYR GQIVQSKINS FRKPLQVKDE SSAATKKLSA TIPKATKPQP VNTSSVTVKS NRSSNMTATT KFVSTTSQNT QLVRPPIRSH HSNTRDTVQK GISRTSANVT IRKGPHEKEL LQSKTALSSV KTSSSQGIIR NKTLRSIAS EVIARPASLS NDKLMEKSEP VDQRRHTAGK AIVDSRSAQP KETSEERKAR LSEWKAGKGR VLKRPPNSVV TQHEPAGQNE KPVGSFWTTM AEDEQRLFT EKVNNTFSEC LNLINEGCPK EDILVTLNDL IKNIPDAKKL VKYWICLALI EPITSPIENI IAIYEKAILA GAQPIEEMRH TIVDILTMKS QEKANLGENM EKSCASKEEV KEVSIEDTGV DVDPEKLEME SKLHRNLLFQ DCEKEQDNKT KDPTHDVKTP NTETRTSCLI KYNVSTTPYL QSVKKKVQFD GTNSAFKELK FLTPVRRSRR LQEKTSKLPD MLKDHYP CVS SLEQLTELGR ETDAFVCRPN AALCRVYYEA DTT
Source	Mammalian cell
Target Names	CKAP2
Protein Names	Recommended name: Cytoskeleton-associated protein 2 Alternative name(s): CTCL tumor antigen se20-10 Tumor- and microtubule-associated protein
Expression Region	1-683
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.