



# Recombinant Human Gamma-crystallin N (CRYGN)

<b>Product Code</b>	CSB-YP841000HU
<b>Abbreviation</b>	CRYGN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WXF5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAQRSGKITL YEGKHFTGQK LEVFGDCDNF QDRGFMNRVN SIHVESGAWV CFNHPDFRGQ QFILEHGDYP DFFRWNSHSD HMGSCRPVGM HGEHFRLEIF EGCNFTGQCL EFLEDSPFLQ SRGWVKNCVN TIKVYGDGAA WSPRSFGAED FQLSSSLQSD QGPEEATTKP ATTQPPFLTA NL
<b>Source</b>	Yeast
<b>Target Names</b>	CRYGN
<b>Protein Names</b>	Recommended name: Gamma-crystallin N Alternative name(s): Gamma-N-crystallin
<b>Expression Region</b>	1-182
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The beta and gamma crystallins are evolutionarily related families of proteins that are localized to the refractive structure of the eye lens. This protein is unique in that it has both beta and gamma crystallin protein motifs.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.