



Recombinant Human THAP domain-containing protein 4 (THAP4)

Product Code	CSB-MP841222HU
Abbreviation	THAP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WY91
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVICCAAVNC SNRQGKGEKR AVSFHRFPLK DSKRLIQWLK AVQRDNWTPT KYSFLCSEHF TKDSFSKRLE DQHRLLKPTA VPSIFHLTEK KRGAGGHGRT RRKDASKATG GVRGHSSAAT SRGAAGWSPS SSGNPMAPKE SRRLKQAALQ GEATPRAAQE AASQEQAAQA LERTPGDGLA TMVAGSQGKA EASATDAGDE SATSSIEGGV TDKSGISMDD FTPPGSGACK FIGSLHSYSF SSKHTRERPS VPREPIDRKR LKKDVEPSCS GSSLGPKDGL AQSPSSSLT ATPQKPSQSP SAPPADVTPK PATEAVQSEH SDASPMISINE VILSASGACK LIDSLHSYCF SSRQNKSQVC CLREQVEKKN GELKSLRQRV SRSDSQVRKL QEKLDELRRV SVPYPSSLLS PSREPPKMNP VVEPLSWMLG TWLSDPPGAG TYPTLQPFQY LEEVHISHVG QPMLNFSFNS FHPDTRKPMH RECGFIRLKP DTNKVAFVSA QNTGVVEVEE GEVNGQELCI ASHSIARISF AKEPHVEQIT RKFRLNSEGK LEQTVSMATT TQPMTQHLHV TYKKVTP
Source	Mammalian cell
Target Names	THAP4
Protein Names	Recommended name: THAP domain-containing protein 4
Expression Region	1-577
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.