



# Recombinant Human THAP domain-containing protein 4 (THAP4)

<b>Product Code</b>	CSB-MP841222HU
<b>Abbreviation</b>	THAP4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WY91
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVICCAAVNC SNRQGKGEKR AVSFHRFPLK DSKRLIQWLK AVQRDNWTPT KYSFLCSEHF TKDSFSKRLE DQHRLLKPTA VPSIFHLTEK KRGAGGHGRT RRKDASKATG GVRGHSSAAT SRGAAGWSPS SSGNPMAPKE SRRLKQAALQ GEATPRAAQE AASQEQAQQA LERTPGDGLA TMVAGSQGKA EASATDAGDE SATSSIEGGV TDKSGISMDD FTTPGSGACK FIGSLHSYSF SSKHTRERPS VPREPIDRKR LKKDVEPSCS GSSLGPDKGL AQSPSSSLT ATPQKPSQSP SAPPADVTPK PATEAVQSEH SDASPMISINE VILSAGACK LIDSLHSYCF SSRQNKSQVC CLREQVEKKN GELKSLRQRV SRSDSQVRKL QEKLDLRRV SVPYPSSLLS PSREPPKMNP VVEPLSWMLG TWLSDPPGAG TYPTLQPFQY LEEVHISHVG QPMLNFSFNS FHPDTRKPMH RECGFIRLKP DTNKVAFVSA QNTGVVEVEE GEVNGQELCI ASHSIARISF AKEPHVEQIT RKFRLNSEGK LEQTVSMATT TQPMTQHLHV TYKKVTP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	THAP4
<b>Protein Names</b>	Recommended name: THAP domain-containing protein 4
<b>Expression Region</b>	1-577
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.