



Recombinant 3-ketoacyl-CoA thiolase (fadA)

Product Code	CSB-MP841272EOD
Abbreviation	fadA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X8J4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MEQVVIVDAI RTPMGRSKGG AFRNVRAEDL SAHLMRSLLA RNPALAAAL DDIYWGCVQQ TLEQGFNIAR NAALLAEVPH SVPVTVNRL CGSSMQALHD AARMIMTGDA QAQLVGGVEH MGHVPM SHGV DFHPGLSRNV AKAAGMMGLT AEMLARMHGI SREMQDAFAA RSHARAWAAT QSGAFKNEII PTGGHDADGV LKQFNIDEVI RPETTVEALA TLRPAFDPVN GTVTAGTSSA LSDGAAAMLV MSESRAHELK LKPRARVRSM AVVGCDPSIM GYGPVPASKL ALKKAGLSVS DIGVFEMNEA FAAQILPCIK DLGLIEQIDE KINLNGGAIA LGHPLGCSGA RISTTLLNLM EHKDVQFGLA TMCIGLGQGI ATVFERV
Source	Mammalian cell
Target Names	fadA
Protein Names	Recommended name: 3-ketoacyl-CoA thiolase EC= 2.3.1.16 Alternative name(s): Acetyl-CoA acyltransferase Beta-ketothiolase Fatty acid oxidation complex subunit beta
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.